

All Purpose Blood Pressure Cuffs

SunTech Medical
Brilliant Blood Pressure Solutions™



All Purpose Cuffs CLINICAL GRADE BP CUFF

Cuffs designed beyond
your expectations

- Designed with the patient in mind
- Customizable
- Latex free
- Meets all US and International regulatory requirements



All Purpose Blood Pressure Cuffs



*Designed beyond
your expectations.*

Ranging from Infant to Thigh, our color-coded All Purpose cuffs are manufactured with easy to clean high quality nylon, which provides an outstanding combination of comfort and durability. These cuffs are easy to use and maintain in the widest variety of clinical applications.

All of our cuffs are designed latex free, feature easily removable bladders to aid in washing and disinfecting, are in compliance with Council Directive 93/42/EEC Medical Devices Directive, and adhere to the recommendations in the AAMI SP10:2002.

FEATURES AND BENEFITS

Ease of Use: Simple application with pictorial indicators, range and artery indicators.











Range of Sizes: Color coded cuffs ranging from Infant to Thigh with pictorial size indicators.


Markets: General BP, patient monitoring.

Customized Labeling: Accommodating special needs by creating partnerships for private labeling with a minimum of 500 cuffs annually.

Cleaning: Sturdy nylon construction makes the cuff easy to clean and disinfect with a removable bladder for machine washing.

Ranges and Symbols

	Infant	8 - 13 cm
	Child	12 - 19 cm
	Child Long	12 - 19 cm
	Small Adult	17 - 25 cm
	Small Adult Long	17 - 25 cm
	Adult	23 - 33 cm
	Adult Long	23 - 33 cm
	Large Adult	31 - 40 cm
	Large Adult Long	31 - 40 cm
	Thigh	38 - 50 cm



Index Line.



Range markers.



Artery symbol.



Circumference.



Product contains no PVC (polyvinyl chloride) material.



Product contains no latex material.



In compliance with Council Directive 93/42/EEC Medical Devices Directive.

