










OPTIONAL REPLACEMENT CONNECTORS			
<i>Photo</i>	<i>Description</i>	<i>SunTech Part Number</i>	<i>Monitors</i>
Screw Type (Threaded)			
	Plastic Male Screw-Type (Threaded), 10/pack	64-0076-11	SunTech 247 GE Critikon Dinamap Criticare 507 Welch Allyn Altas Welch Allyn Vital Signs Welch Allyn Spot Spacelabs Nihon Kohden (Lifescope L)
Luer Slips			
	Metal Male Luer Slip, 10/pack	64-0047-15	
	Plastic Male Luer Slip, 10/pack	64-0123-10	CASMED 740/750
	Plastic Female Luer Slip, 10/pack	64-0125-10	Baum
Quick Connect			
	Plastic Male Quick Connect, 10/pack	64-0047-13	Nihon Kohden Datex-Ohmeda*
	Plastic Female Quick Connect, 10/pack	64-0047-12	SunTech Tango+ SunTech Oscar 2 GE/Marquette two-tube Datex-Ohmeda* Protocol Propaq Welch Allyn Propaq
Bayonet (Rectus)			
	Metal Bayonet (Rectus), 10/pack	64-0047-14	Philips/HP Siemens Colin Spacelabs Datascope (Accurtor Plus & Passport 2)
	Plastic 10/pack	64-0122-10	Philips/HP Siemens Colin Spacelabs Datascope (Accurtor Plus & Passport 2)
Inline			
	Plastic Female Inline, 10/pack	64-0124-10	Marquette

*Datex-Ohmeda uses one male and one female quick connect on its two-tube monitors