

# Eclipse™ Blood Pressure Cuffs



Cuffs designed beyond  
your expectations

- Designed with the patient in mind
- Customizable
- Latex free
- Meets all US and International regulatory requirements



# Eclipse™ Blood Pressure Cuffs



*Designed beyond  
your expectations.*

Designed with unique features that simplify cuff position and reduce the possibility of improper use, our Eclipse D-Ring cuffs are a great choice for home monitoring, all self-application environments and anywhere that speed and accuracy are important. Our Adult and Large Adult sizes are curved for a more tailored fit.

All of our cuffs are designed latex free, feature easily removable bladders to aid in washing and disinfecting, are in compliance with Council Directive 93/42/EEC Medical Devices Directive, and adhere to the recommendations in the AAMI SP10:2002.

## FEATURES AND BENEFITS

<b>Patient Comfort:</b>	Innovative design that prevents pinching and brings greater patient comfort to blood pressure monitoring.	<b>Latex Free:</b>	Worry free use on all patients.
<b>Range of Sizes:</b>	A full range of sizes to accommodate Pediatric through Large Adult. The Adult and Large Adult sizes are curved for a more tailored fit.	<b>Markets:</b>	Home BP units, self-application environments, general BP.
<b>Ease of Application:</b>	D-Ring design makes application easier for the practitioner and patient.	<b>Customized Labeling:</b>	Accommodating special needs by creating partnerships for private labeling with a minimum of 500 cuffs annually.
		<b>Cleaning:</b>	Sturdy nylon construction makes the cuff easy to clean and disinfect with a removable bladder for hand washing.

## Ranges and Symbols

Pediatric	16 - 22 cm
Small Adult	21 - 29 cm
Adult	28 - 37 cm
Large Adult	36 - 46 cm



Range markers.

ARTERY



Artery symbol.



Circumference.



Product contains no PVC (polyvinyl chloride) material.



Product contains no latex material.