

Orbit™ Blood Pressure Cuffs

 SunTech Medical®
Brilliant Blood Pressure Solutions™



Orbit™
BP CUFF

Cuffs designed beyond
your expectations

- Designed with the patient in mind
- Ease of application
- Latex free
- Meets all US and international regulatory requirements



Orbit™ Blood Pressure Cuffs



*Designed beyond
your expectations.*













Slip on an Orbit cuff and you will appreciate how comfortable a blood pressure cuff can be. More than any other product on the market, the Orbit's state of the art design adds comfort for patients who must wear a cuff for an extended period of time. The soft, supple sleeve, made of lycra/nylon stretch material, gently hugs the arm holding the cuff in place. Our Adult Plus and Large Adult sizes are curved for a more tailored fit. All this technology provides uninterrupted accurate measurements even when patients are on the move.

All of our cuffs are designed latex free, feature easily removable bladders to aid in washing and disinfecting, are in compliance with Council Directive 93/42/EEC Medical Devices Directive, and adhere to the recommendations in the AAMI SP10:2002.

FEATURES AND BENEFITS

Patient Comfort:	Innovative sleeve design that secures cuff to the arm and provides comfort for longer testing periods.	Latex Free:	Worry free use on all patients.
Range of Sizes:	A full range of sizes from Small Adult to Large Adult. The Adult and Large Adult sizes are curved for a more tailored fit.	Markets:	Ambulatory BP, extended-wear testing, self-application environments.
Ease of Application:	Pictorial application instructions, artery markers and clear range indicators make it easy to apply the cuff.	Cleaning:	Sturdy material and a removable bladder for machine washing.

Ranges and Symbols

<table border="0"> <tr> <td>Small Adult</td> <td>18 - 27 cm</td> </tr> <tr> <td>Adult</td> <td>25 - 35 cm</td> </tr> <tr> <td>Adult Plus</td> <td>33 - 40 cm</td> </tr> <tr> <td>Large Adult</td> <td>39 - 46 cm</td> </tr> </table>	Small Adult	18 - 27 cm	Adult	25 - 35 cm	Adult Plus	33 - 40 cm	Large Adult	39 - 46 cm	  	<table border="0"> <tr> <td>Range markers.</td> <td></td> <td>Product contains no PVC (polyvinyl chloride) material.</td> </tr> <tr> <td>Artery symbol.</td> <td></td> <td>Product contains no latex material.</td> </tr> <tr> <td>Circumference.</td> <td></td> <td>In compliance with Council Directive 93/42/EEC Medical Devices Directive.</td> </tr> </table>	Range markers.		Product contains no PVC (polyvinyl chloride) material.	Artery symbol.		Product contains no latex material.	Circumference.		In compliance with Council Directive 93/42/EEC Medical Devices Directive.
Small Adult	18 - 27 cm																		
Adult	25 - 35 cm																		
Adult Plus	33 - 40 cm																		
Large Adult	39 - 46 cm																		
Range markers.		Product contains no PVC (polyvinyl chloride) material.																	
Artery symbol.		Product contains no latex material.																	
Circumference.		In compliance with Council Directive 93/42/EEC Medical Devices Directive.																	